



PennState

**One Health
Microbiome Center**

KickStart Workshop

Fall 2026 KickStart Workshop for the analysis of microbial communities

Day 2: Amplicon Data Analysis with Qiime2

Veronica Roman-Reyna, PhD

Assistant Professor

Plant Pathology and Environmental Microbiology

roman-reyna@psu.edu

Day 2

- About amplicon sequencing
- Data analysis with Qiime
- Lunch break
- Cont. data analysis with Qiime

Learning Goals



Concepts

What is an amplicon, why is a good barcode, and how sequencing it reshaped our view of microbial diversity.



Tools

Which software options exist for amplicon analysis, and how to work with QIIME 2.



Data

Data and databases. NCBI SRA.




Workflow

The QIIME 2 workflow: from raw reads to diversity, taxonomy, and differential-abundance results.

What is an amplicon?

- A DNA fragment produced by PCR **amplification** of one specific, **targeted** region of a genome.
- In microbiome work, that target is a taxonomically informative marker gene.
- Millions of copies of that one region, from every organism in the sample, are amplified and sequenced together.




Taxonomically informative marker:

The most common markers are Ribosomal RNA (rRNA)

What is the S?
Svedberg unit.

A measure of how fast a particle sediments in a centrifuge.




1990 vs 2016

Evolution: Woese *et al.*

16S*


Proc. Natl. Acad. Sci. USA 87 (1990)




Sources:

<https://doi.org/10.1073/pnas.87.12.4576>


<https://doi.org/10.1038/nmicrobiol.2016.48>




16s rRNA structure



Nature Reviews | Microbiology



Slide modified from J. Wan




Amplicon sequencing

OTUs	ASVs or ESVs
Operational Taxonomic Unit	Amplicon/Exact Sequence Variant
The main differences are resolution, methodology, reproducibility.	
1 OTU $\geq$ 97%	1 ASV = single nucleotide differences.

There has been a push since 2017 to only use ASVs (<https://www.nature.com/articles/ismej2017119>)

Amplicon sequencing



Amplicon Sequencing: Strengths & Limitations

When to use it

- Questions about diversity and/or compositional shifts.
- You need to process many samples cost-effectively (~\$20* vs ~\$150*)
- Low abundance organisms (PCR)

Bacteria and Archaea (16S rRNA)

Fungi (ITS, 18S rRNA)

Eukaryotes (18S rRNA)

When not to

- Species/Strain-level resolution or functional gene content.
- Gene copy number varies across taxa.

DNA/RNA virus

Archaea, bacteria, Eukaryotes

Software for Amplicon Analysis



QIIME 2

A plugin-based platform covering the whole workflow (Conda and terminal).



DADA2 (R package)

Directly in R for users who want full scripting control outside QIIME 2.



Mothur

One of the original open-source amplicon pipelines. **OTU**-based workflows.



- Quantitative insights into Microbial Ecology
- Pronounced like "chime".
- Created ~2010 during the Human Microbiome Project (2007 - 2016) under the leadership of Greg Caporaso and Rob Knight.
- Essentially, QIIME is a set of **commands** to transform microbiome **data** into **intermediate outputs** and **visualizations**.



QIIME 1 is outdated and QIIME 2 contains many improvements to analytic functionality, is easier to install, and is easier to use.

[nature](#) > [nature methods](#) > [correspondence](#) > [article](#)

Correspondence | [Published: 11 April 2010](#)

QIIME allows analysis of high-throughput community sequencing data

[J Gregory Caporaso](#), [Justin Kuczynski](#), [Jesse Stombaugh](#), [Kyle Bittinger](#), [Frederic D Bushman](#), [Elizabeth K Costello](#), [Noah Fierer](#), [Antonio Gonzalez Peña](#), [Julia K Goodrich](#), [Jeffrey I Gordon](#), [Gavin A Huttley](#), [Scott T Kelley](#), [Dan Knights](#), [Jeremy E Koenig](#), [Ruth E Ley](#), [Catherine A Lozupone](#), [Daniel McDonald](#), [Brian D Muegge](#), [Meg Pirrung](#), [Jens Reeder](#), [Joel R Sevinsky](#), [Peter J Turnbaugh](#), [William A Walters](#), [Jeremy Widmann](#), ... [Rob Knight](#)  [+ Show authors](#)

[Nature Methods](#) **7**, 335–336 (2010) | [Cite this article](#)

67k Accesses | **25k** Citations | **584** Altmetric | [Metrics](#)

[nature](#) > [nature biotechnology](#) > [correspondence](#) > [article](#)

Correspondence | [Published: 24 July 2019](#)

Reproducible, interactive, scalable and extensible microbiome data science using QIIME 2

[Evan Bolyen](#), [Jai Ram Rideout](#), [Matthew R. Dillon](#), [Nicholas A. Bokulich](#), [Christian C. Abnet](#), [Gabriel A. Al-Ghalith](#), [Harriet Alexander](#), [Eric J. Alm](#), [Manimozhiyan Arumugam](#), [Francesco Asnicar](#), [Yang Bai](#), [Jordan E. Bisanz](#), [Kyle Bittinger](#), [Asker Brejnrod](#), [Colin J. Brislawn](#), [C. Titus Brown](#), [Benjamin J. Callahan](#), [Andrés Mauricio Caraballo-Rodríguez](#), [John Chase](#), [Emily K. Cope](#), [Ricardo Da Silva](#), [Christian Diener](#), [Pieter C. Dorrestein](#), [Gavin M. Douglas](#), ... [J. Gregory Caporaso](#)  [+ Show authors](#)

[Nature Biotechnology](#) **37**, 852–857 (2019) | [Cite this article](#)

77k Accesses | **8665** Citations | **252** Altmetric | [Metrics](#)

Amplicon

A suite of plugins that provide broad analytic functionality to support microbiome marker gene analysis from raw sequencing data through publication quality visualizations and statistics.

Version

2026.1

- Links
- [Documentation](#)[Source Code](#)[Conda Channel](#)

Built-in plugins

[q2](#), [q2-alignment](#), [q2-boots](#), [q2-composition](#), [q2-cutadapt](#), [q2-dada2](#), [q2-deblur](#), [q2-demux](#), [q2-diversity](#), [q2-diversity-lib](#), [q2-emperor](#), [q2-feature-classifier](#), [q2-feature-table](#), [q2-fondue](#), [q2-fragment-insertion](#), [q2-kmerizer](#), [q2-longitudinal](#), [q2-metadata](#), [q2-mystery-stew](#), [q2-phylogeny](#), [q2-quality-control](#), [q2-quality-filter](#), [q2-sample-classifier](#), [q2-stats](#), [q2-taxa](#), [q2-types](#), [q2-vizard](#), [q2-vsearch](#), [rachis](#), [RESCRIPT](#)

Qiime offers multiple tutorials and tools explanations

qiime2docs

Search

K

Microbiome marker gene analysis with QIIME 2

Tutorials>How To Guides>Explanations>References>Back Matter>

Microbiome marker gene analysis with QIIME 2

Welcome! 🙌 This is the primary documentation introducing the use of QIIME 2 for marker gene (i.e., amplicon) based microbiome analysis.

Transition from “the old docs”

As of April 2025, this site replaces the old QIIME 2 user documentation, <https://docs.qiime2.org>. We’ve prioritized content to transition from the old documentation based on our website analytics, so the most frequently accessed content is already here. If you’re looking for content from the old QIIME 2 user documentation, you can find it at <https://docs.qiime2.org/2024.10/> (but please consider [letting us know](#) what that content is, as we’re trying to transition everything that’s needed to this site).

CONTENTS

Navigating this documentationContributorsFunding 🙌Citing QIIME 2License

Qiime2 - What do you need



Using conda for installation.

“analogy”

conda installations



source code*



Sequencing reads for analysis

Type of reads

fastq/fq files → Raw reads. Files that contains nucleotide and their quality.
Single-end → One file per sample (sample.fq)
[Seq. technologies: Illumina, PacBio, Oxford Nanopore]
Paired-end → A pair of files per sample (sample_R1.fq, sample_R2.fq).
[Seq. technologies: Illumina]

Source

- Your own data
- Publicly available data from tutorials**
- The Sequence Read Archive stores raw reads from published studies.

Once you have ASVs, you need a reference to assign taxonomy.

SILVA

16S/18S rRNA database spanning Bacteria, Archaea, and Eukaryotes. Updated in 2024 (<https://www.arb-silva.de/>).

GTDB

Genome Taxonomy Database, taxonomy standardized by genome phylogeny (<https://gtdb.ecogenomic.org/>).

***Greengenes2**

A 16S rRNA database re-built on a genome-based reference tree. Updated in 2024.

UNITE


The standard reference for fungal ITS sequences. (<https://unite.ut.ee/>).

NCBI

The largest general sequence archive. You can build a custom reference from it (<https://www.ncbi.nlm.nih.gov/>).

- QIIME 2 pre-trained classifiers for several of these are available at Qiime website.
- Check which version of Qiime you are using as that will affect what database you download.

The NCBI SRA (Sequence Read Archive)



What is SRA?

“The [Sequence Read Archive](#) (SRA) is the [National Center for Biotechnology Information](#) (NCBI) database that stores sequence data obtained from next generation sequence technology.”

How is SRA useful?

Deposit your own data while submitting an article, Data mining (use data to ask different questions)

How to use SRA?

Online and command line (SRA tool kit: <https://github.com/ncbi/sra-tools>)

Full

Send to

SRX22919344: MySeq of soil microbiome

1 ILLUMINA (Illumina MiSeq) run: 116,635 spots, 69.6M bases, 40.8Mb downloads

Design: The metagenomic DNA from each soil sample was extracted using the FastDNA spin Kit for soil from MP biomedical, following the manufacturers instruction. DNA concentration was determined by Nanodrop (Thermo Fisher) and the V3-V4 hypervariable 16s rRNA region (primers Pro341F 5-CCTACGGGNNBGCASCAG -3 and Pro805R 5-GACTACNVGGGTATCTAATCC-3) was amplified. 16S rRNA library preparation was performed by BMR Genomics (Italy) for Illumina MiSeq 300PE sequencing.

Submitted by: University of Perugia

Study: Soil microbiota of organic and conventional management

[PRJNA1053916](#) • [SRP478556](#) • [All experiments](#) • [All runs](#)
[show Abstract](#)

Sample: Conventional low-input

[SAMN38879569](#) • [SRS19887352](#) • [All experiments](#) • [All runs](#)
Organism: [soil metagenome](#)

Library:

Name: 5A_Low_R2_G
Instrument: Illumina MiSeq
Strategy: AMPLICON
Source: METAGENOMIC
Selection: PCR
Layout: PAIRED

Runs: 1 run, 116,635 spots, 69.6M bases, [40.8Mb](#)

Run	# of Spots	# of Bases	Size	Published
SRR27240947	116,635	69.6M	40.8Mb	2023-12-18

Recent activity

[Turn Off](#) [Clear](#)

 [soil metagenome \(818729\)](#)

SRA

 [metagenome \(4013063\)](#)

SRA

[See more...](#)

Full ▾

SRX22885110: Soil metagenome

1 DNBSEQ (DNBSEQ-T7) run: 33.5M spots, 10G bases, 6Gb downloads

Design: normal metagenome of soil**Submitted by:** Sichuan Normal University**Study:** Soil metagenome Raw sequence reads[PRJNA1051770](#) • [SRP477915](#) • [All experiments](#) • [All runs](#)[show Abstract](#)**Sample:** T3[SAMN38796427](#) • [SRS19857384](#) • [All experiments](#) • [All runs](#)*Organism:* [soil metagenome](#)**Library:***Name:* S427*Instrument:* DNBSEQ-T7*Strategy:* WGS*Source:* METAGENOMIC*Selection:* RANDOM*Layout:* PAIRED**Runs:** 1 run, 33.5M spots, 10G bases, [6Gb](#)

Run	# of Spots	# of Bases	Size	Published
SRR27205348	33,484,832	10G	6Gb	2023-12-14


ID: 30942389

Check point

- What is an amplicon and when to use it
- Why we use ribosomal sequences as markers
- What is OTU and ASV
- What software can we use
- Work with Qiime2



The QIIME 2 Pipeline



The QIIME 2 Pipeline – Go to the terminal

Home Directory

Work

/storage/group/one/default

/storage/group/one/default

/storage/group/srb6251/default

↑ / storage / home / vfr5124 /

Change directory

Copy path

☒ Show Owner/Mode ☐ Show Dotfiles Filter:

Showing 30 of 65 rows - 0 rows selected

☐	Type ▲	Name		Size	Modified at	Owner	Mode
☐	📁	bin	⋮	-	2/25/2026 10:15:16 AM	vfr5124	775
☐	📁	Desktop	⋮	-	8/5/2025 3:22:56 PM	vfr5124	755

```
***** WARNING *****
This computer system is the property of the Pennsylvania State University. It
is for authorized use only. By using this system, all users acknowledge notice
of, and agree to comply with, the Institute for Computational and Data Science
Account Usage Policies, Data Policies, and relevant University Policies.

Unauthorized or improper use of this system may result in administrative
disciplinary action, civil charges/criminal penalties, and/or other sanctions
as set forth in University Policies. By continuing to use this system
you indicate your awareness of and consent to these terms and conditions of
use.

LOG OFF IMMEDIATELY if you do not agree to the conditions stated in
this warning.

NOTICE: Microsoft MFA enabled users will get a silent push to your two factor
device. After you enter your password please check your two factor device.

***** WARNING *****
[vfr5124@submit02 ~]$
```